

ENVIRONICS TECHNICAL BULLETIN

SERIES 2000/9100 RS232 CABLE CONNECTIONS

The Series 2000/9100 systems include two RS232 male DB25 connectors.
A null modem cable wired like shown below is required for connection to a PC.

| Signal | PC (DB9) | PC (DB25) | | Series 2000 (DB25) | Signal |
|---------------|-----------------|------------------|--|---------------------------|---------------|
| Receive (RD) | 2 | 3 | | 2 | Transmit (TD) |
| Transmit (TD) | 3 | 2 | | 3 | Receive (RD) |
| Ground (GND) | 5 | 7 | | 7 | Ground (GND) |
| | | | | | |

Other signals are not required.