

MODEL 6100 QUICK START PROCEDURE

This procedure outlines the basic steps to setup a blend with or without ozone.

1. Connect power and gases to the system as prescribed in the 6100. See Table of Contents section, "Power and Gas Connections" in the User's Guide and Technical manual.
2. Enter information in GASES mode. Identify, to the system, the gas types and concentrations connected to the rear panel ports of the system. Reference "Gas Setup Mode" in 6100 User's manual.
3. Setup the blend routine in the CONC mode. Reference the "Concentration Mode" in the 6100 User's manual.