

Installing the Environics S4000 software with Windows 7

*Note - This applies to Windows 7 - 32 bit only.
Due to limitations of XP Mode with Windows 7 - 64 bit,
the S4000 software cannot be installed.*

If you are running Windows 7, you will need to download and install "XP Mode" first, before installing the S4000 software.

Click **Start - Windows Virtual PC - Windows XP Mode**

*If Windows XP Mode is not listed in the Start menu, it has not been pre-installed.
You can download it for free from Microsoft.com
<http://www.microsoft.com/windows/virtual-pc/download.aspx>*

Insert the Environics CD

Inside the XP Mode window, click the Start button

Select My Computer, double click the CD or DVD icon.

Locate the SETUP.EXE file on the DVD and double click it to start installation of the S4000 software

After software is installed, go to the top of the XP Mode windows, select Tools - Settings.

Make sure COM1 is set to the Physical COM port on the PC.

If you have any questions, please feel free to contact us at (860) 872-1111 or online at <http://www.environics.com/service-request/>.
